Introduction to Cell Biology Section 1  Lab Quiz #2  Bacterial Staining
1. Give two staining methods used for Petri dishes. (4 points)
a. Comassie flue
b. Gel Code Dhur
2. Give three staining methods used for slides. (5 points)
a. Methyleneblue
b. Crystal/Gram Stain
c. Manthair India Ink
3. What two stains are used in conjunction to distinguish between gram positive and gram negative bacteria? (2 points)
methylexeblue and crystal violet
Methylere blue and crystal violet  3. E.coli are Gram-Maative and will stain while L.P are Gram- positive and will be colored. (4 points)
5. Coomassie Blue forms bonds with NH3+ groups. Why is this stain important? (5 points)
6. Explain how and why we heat fix bacteria. (5 points) wateria is staining
so that the stains don't wash the bacteria

Off the Slide-